INSTRUCTIONS FOR USE
LACRIFILL® Canalicular Gel
Disposable syringe
Sterile Canalicular gel 20mg/ml

I. DESCRIPTION & COMPOSITION

LacriFILL® is sterile, colorless gel composed of non-animal hyaluronic acid (HA) cross-linked with 1,4- butanedioldiglycidylether (BDDE) and reconstituted in a physiological buffer. The product is supplied in 1mL syringe (filling volume: 0.6mL), sterilized using moist heat sterilization process.

The pre-filled syringe is assembled with a plunger rod and a finger rest to allow gel placement. It is then packed in a secondary packaging (blister) to protect it from damage during transportation and storage and corresponds to one unit.

LacriFILL® is available in boxes of 1 or 10 units. **LacriFILL®** does not include any accessory. Qualitative and quantitative components of the injected gel:

Ingredients	Content per g
Crosslinked sodium hyaluronate	20 mg
Sodium chloride	9 mg
Potassium dihydrogen phosphate, KH ₂ PO ₄	0.02 mg
Disodium hydrogen phosphate dihydrate, Na ₂ HPO ₄ · 2H ₂ O	0.22 mg
Water for injection (WFI)	ad 1 g

II. INTENDED USE

a/ Intended purpose and indication

LacriFILL® is designed to temporarily block tear drainage by the occlusion of the lacrimal canaliculus of one or both eyelids in a given patient, thus maintaining lubricating tears on the surface of the eye.

LacriFILL® is intended for use, for up to 6 months, in adult patients experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations, which can be relieved by blocking of the canicular system.

b/ Intended users

LacriFILL® shall only be administered by appropriately trained healthcare professionals who are qualified or accredited in accordance with national law, typically ophthalmologists or optometrists.

c/ Target population

LacriFILL® is used in adults experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations.

d/ Use environment

LacriFILL® is intended to be used in clinics/ hospitals or private practices.

III. CONTRAINDICATIONS

LacriFILL® is contraindicated for patients with lacrimal outflow obstruction, inflammation of the ocular surface not associated with dry eye, active ocular or periorcular infection (for example, dacryocystitis, canaliculitis, conjunctivitis, keratitis) and those who are allergic to hyaluronic acid.

IV. MODE OF ACTION

LacriFILL® is administered into the lower canaliculus of the patient's eyelid using a commercially available LacriFILL canalicula attached to the 0.6mL pre-filled **LacriFILL®** syringe. **LacriFILL®** is designed to temporarily block (for up to 6 months) tear drainage by the occlusion of the lacrimal canalicula in a given patient, thus maintaining lubricating tears on the surface of the eye.

V. PREREQUISITES BEFORE USE AND OPERATING INSTRUCTIONS

Use of sterile goggles is recommended.

Canula Preparation

- A sterile, disposable or reusable 25–27-gauge lacrimal cannicula with Luer-lock connection may be used with **LacriFILL®** injection device. If a reusable lacrimal cannicula, is to be used, the cannicula should be cleaned and sterilized prior to use for insertion of the gel containing with the cannicula's instructions for use and institutional procedures. Depending on the use of 25 or 27-gauge lacrimal cannicula, the extrusion force will be slightly different, requiring more force to provide the same amount of gel using a 27-gauge lacrimal cannicula.
- Prior to punctal irrigation and insertion of **LacriFILL®**, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort. Punctal irrigation must demonstrate good lacrimal drainage prior to insertion of **LacriFILL®**.
- To prepare for insertion of the gel, open the box, inspect the packaging for any signs of damage, remove the tap containing the gel-filled syringe, and peel off tap cover. Do not use if packaging has been compromised or device has been damaged.
- Remove the cap from the gel-filled syringe and attach a sterile lacrimal cannicula to the syringe (Figure 1).

Insertion

- Place the tip of the lacrimal cannicula into the lower punctum of the first eyelid (Figure 4) and insert 0.2 mL of gel (Figure 5). The punctum may be dilated if necessary for introduction of the lacrimal cannicula. Note: In approximately one-third of cases, gel extrusion from the upper punctum might be observed in which case excess gel can be irrigated from the ocular surface if it obscures the patient's vision.


Figure 4: Place cannicula tip in lower punctum


Figure 5: Deliver 0.2 mL gel into lower canaliculus

- For insertion of the gel into the punctum of the second eye, purge the cannicula with 0.1 mL of gel (the volume of gel left in the cannicula should be 0.2 mL) then repeat as illustrated in Figures 4 and 5 above.
- Remove the lacrimal cannicula from the syringe. If using a reusable cannicula, clean and re-sterilize using standard procedures consistent with the cannicula manufacturer's instructions for use and institutional procedures.
- Discard the syringe, the disposable cannicula if used, any remaining gel, and packaging in accordance with accepted medical practice and applicable local, and national requirements.
- Used the syringe, the disposable cannicula if used, any remaining gel, and packaging in accordance with accepted medical practice and applicable local, and national requirements.

Removal of the Gel

In case of any adverse event, **LacriFILL®** may be removed. To remove the gel, irrigate the gel from the cannicula using a syringe filled with sterile saline. Attach a lacrimal cannicula onto the syringe, insert the cannicula into the canaliculus, and gently irrigate until **LacriFILL®** is removed from the canaliculus as evidenced by free flow of fluid through the lacrimal drainage system and/or the patient's report of swallowing or a fluid sensation in the nose or throat.

VI. WARNINGS

- As with any procedure, use of **LacriFILL®** carries a risk of infection. Standard precautions associated with injectable materials should be followed. Should an infection occur, flush the gel from the canaliculus.
- If the patient experiences irritation, infection, or epiphora after insertion of **LacriFILL®**, the gel should be removed.
- In the pivotal clinical study, the safety and effectiveness of **LacriFILL®** was characterized over a 6-month period. The safety and effectiveness of the device for longer periods of use has not been established. Therefore **LacriFILL®** should be used for no more than 6 months. After 6 months, **LacriFILL®** should be removed by thorough lacrimal irrigation.
- When removing **LacriFILL®**, ensure that punctal irrigation is sufficiently thorough. If irrigation does not remove the gel, lacrimal probing or surgical exploration may be needed.

- To minimize the risks of potential complications, **LacriFILL®** should only be inserted by health care practitioners who have appropriate training and experience, and who are knowledgeable about the anatomy at and around the punctum and canaliculus.
- Obstruction of the canaliculus should be determined prior to insertion of **LacriFILL®**.
- Do not insert **LacriFILL®** if the canaliculus is obstructed. Care should be used not to perforate the canaliculus during the insertion of the **LacriFILL®**. Perforation may cause pain and increased risk of infection. If perforation occurs, delay the insertion of the gel until the wound heals.

VII. PRECAUTIONS

LacriFILL® is a sterile, colorless gel without visible particulate. If the contents of a syringe show signs of separation and/or appear cloudy or colored, do not use the syringe; notify **Nordic Pharma** immediately at lacriFILL-vigilance@nordicpharma.com.

- LacriFILL®** is provided STERILE for single patient use only, for insertion in one or two lacrimal canaliculi during same session to avoid risk of infection.
- Do not re-sterilize the syringe containing the gel.
- Do not insert more or less than 0.2 mL of gel into each lower punctum unless the canalicular anatomy of a given patient requires it.
- Use the device prior to the expiration date printed on the label.
- Prior to use, inspect the package and device to verify that there is no damage to the device or package. Do not use the device if the device or packaging appears open or damaged.
- Failure to comply with the cannicula attachment instructions could result in cannicula disengagement and/or product leakage at the Luer lock and cannicula hub connection.
- Prior to punctal irrigation and insertion of **LacriFILL®**, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort.
- Do not prime more or less than 0.1 mL of gel before insertion in each lower punctum.
- Inadequate or excessive ejection force, or improper technique used to insert **LacriFILL®** into the lower canaliculus may cause patient discomfort.
- After use, dispose of the syringe in accordance with accepted medical practice and applicable local, and state requirements.
- If a reusable lacrimal cannicula has been used, it should be cleaned and re-sterilized using standard procedures consistent with the cannicula manufacturer's instructions for use and institutional procedures.
- As with any punctal or canalicular plug, **LacriFILL®** may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be adjusted accordingly.
- The clinical performance of **LacriFILL®** has not been established in the following conditions/ patient populations:
 - Pregnancy and breastfeeding.
 - Pediatric patients.
 - Use of ophthalmic cyclosporine (Restasis) within 6 months or liftegrast (Xidra) within 3 months.
 - Corneal transplant.
 - Ocular surgery (such as cataract surgery or LASIK) within six months.
- Current systemic infection, uncontrolled ocular surface disease, uncontrolled immunodeficiency disease.
- Significant corneal or conjunctival scarring, pterygium or nodular pinguecula; current ocular infection (except mild blepharitis), conjunctivitis or inflammation not associated with dry eye; anterior (epithelial) basement membrane corneal dystrophy or other clinically significant corneal dystrophy or degeneration; history of ocular herpetic infection; evidence of keratoconus; lid or lacrimal cancer.
- Active severe systemic allergy, seasonal allergies, rhinitis or sinusitis requiring treatment (i.e., antihistamines, decongestants, oral or aerosol steroids).
- Use of steroids, including administration by systemic or topical ocular routes.
- Ectropion.

X. EXPIRATION DATE, AND STORAGE, HANDLING AND DISPOSAL

LacriFILL® has a shelf life of 24 months. **LacriFILL®** must be used prior to the expiration date on the package. Store at a temperature of +5°C up to +25°C. Protect from sunlight. Refrigeration is not required.

After opening, discard any unused material of **LacriFILL®** according to the institutional policies. Do not re-sterilize the syringe containing the gel.

Do not insert more or less than 0.2 mL of gel into each lower punctum unless the canalicular anatomy of a given patient requires it.

Use the device prior to the expiration date printed on the label.

Prior to use, inspect the package and device to verify that there is no damage to the device or package. Do not use the device if the device or packaging appears open or damaged.

Failure to comply with the cannicula attachment instructions could result in cannicula disengagement and/or product leakage at the Luer lock and cannicula hub connection.

Prior to punctal irrigation and insertion of **LacriFILL®**, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort.

Do not prime more or less than 0.1 mL of gel before insertion in each lower punctum.

Inadequate or excessive ejection force, or improper technique used to insert **LacriFILL®** into the lower canaliculus may cause patient discomfort.

After use, dispose of the syringe in accordance with accepted medical practice and applicable local, and state requirements.

If a reusable lacrimal cannicula has been used, it should be cleaned and re-sterilized using standard procedures consistent with the cannicula manufacturer's instructions for use and institutional procedures.

As with any punctal or canalicular plug, **LacriFILL®** may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be adjusted accordingly.

The clinical performance of **LacriFILL®** has not been established in the following conditions/ patient populations:

- Itching
- Local infection
- Local irritation
- Loss / Partial extrusion
- Punctal granuloma
- Swelling

IX. HOW LACRIFILL IS SUPPLIED AND STERILE CONDITIONS

LacriFILL® is supplied in a syringe containing 0.6 mL of gel for possible insertion of both eyes of a single patient during the same session. The contents of the syringe are sterile by moist heat method. Do not re-sterilize. Do not use if packaging is open or damaged.

X. EXPIRATION DATE, AND STORAGE, HANDLING AND DISPOSAL

LacriFILL® has a shelf life of 24 months. **LacriFILL®** must be used prior to the expiration date on the package. Store at a temperature of +5°C up to +25°C. Protect from sunlight. Refrigeration is not required.

After opening, discard any unused material of **LacriFILL®** according to the institutional policies. Do not re-sterilize the syringe containing the gel.

Do not insert more or less than 0.2 mL of gel into each lower punctum unless the canalicular anatomy of a given patient requires it.

Use the device prior to the expiration date printed on the label.

Prior to use, inspect the package and device to verify that there is no damage to the device or package. Do not use the device if the device or packaging appears open or damaged.

Failure to comply with the cannicula attachment instructions could result in cannicula disengagement and/or product leakage at the Luer lock and cannicula hub connection.

Prior to punctal irrigation and insertion of **LacriFILL®**, a topical anesthetic agent should be instilled in the conjunctival sac for patient comfort.

Do not prime more or less than 0.1 mL of gel before insertion in each lower punctum.

Inadequate or excessive ejection force, or improper technique used to insert **LacriFILL®** into the lower canaliculus may cause patient discomfort.

After use, dispose of the syringe in accordance with accepted medical practice and applicable local, and state requirements.

If a reusable lacrimal cannicula has been used, it should be cleaned and re-sterilized using standard procedures consistent with the cannicula manufacturer's instructions for use and institutional procedures.

As with any punctal or canalicular plug, **LacriFILL®** may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be adjusted accordingly.

The clinical performance of **LacriFILL®** has not been established in the following conditions/ patient populations:

As with any procedure, use of **LacriFILL®** carries a risk of infection. Standard precautions associated with injectable materials should be followed. Should an infection occur, flush the gel from the canaliculus.

If the patient experiences irritation, infection, or epiphora after insertion of **LacriFILL®**, the gel should be removed.

In the pivotal clinical study, the safety and effectiveness of **LacriFILL®** was characterized over a 6-month period. The safety and effectiveness of the device for longer periods of use has not been established. Therefore **LacriFILL®** should be used for no more than 6 months. After 6 months, **LacriFILL®** should be removed by thorough lacrimal irrigation.

When removing **LacriFILL®**, ensure that punctal irrigation is sufficiently thorough. If irrigation does not remove the gel, lacrimal probing or surgical exploration may be needed.

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Nordic Group B.V.
Symatose SAS